



Product Datasheet

alpha-2 antiplasmin antibody (orb73687)


Description	alpha-2 antiplasmin antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P08697
Expiration Date	12 months from date of receipt.




Western Blot analysis of A: Hela whole ce...



Immunohistochemical analysis of formalin...



Costained with Hoechst 33342



Confocal immunofluorescence analysis of ...