



Product Datasheet

ELP3 antibody (orb73723)


Description	ELP3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9H9T3
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:NT2D1, B:IMR3...



Merged with DNA probe



Confocal immunofluorescence analysis of ...