

Product Datasheet

ABF3 antibody (orb737390)

Description	Rabbit polyclonal antibody to ABF3
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	Unconjugated
Tested Applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9M7Q3
Dilution Range	WB: 1/500 - 1/2000
Expiration Date	12 months from date of receipt.