


## Product Datasheet

### Anti-ABI5 Antibody (orb737392)

|                            |                                                                                                                                        |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to ABI5                                                                                                     |
| <b>Species/Host</b>        | Rabbit                                                                                                                                 |
| <b>Reactivity</b>          | A. thaliana                                                                                                                            |
| <b>Conjugation</b>         | Unconjugated                                                                                                                           |
| <b>Tested Applications</b> | WB                                                                                                                                     |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| <b>Note</b>                | For research use only                                                                                                                  |
| <b>Clonality</b>           | Polyclonal                                                                                                                             |
| <b>Clone Number</b>        | ABI5                                                                                                                                   |
| <b>Uniprot ID</b>          | <b>Q9SJN0</b>                                                                                                                          |
| <b>Dilution Range</b>      | WB: 1/500 - 1/2000                                                                                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                        |



Western blot analysis of ABI5 expression in Arabidopsis thaliana (A) whole cell lysates.

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)