



Product Datasheet

TAP antibody (orb73740)


Description	TAP antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9UBU9
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:Hela, B:Molt-...



Immunohistochemical analysis of formalin...



Confocal immunofluorescence analysis of ...