

Product Datasheet

Human Serpin I2 protein (orb737423)

Description	Recombinant Human Serine Protease Inhibitor-clade I2 is produced by our Mammalian expression system ...
Endotoxins	Less than 0.1 EU/ug
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	C-6His
Note	For research use only
Uniprot ID	075830
Expiration Date	6 months from date of receipt.