

Product Datasheet

Human Mucin-1 protein (orb737435)

| | |
|------------------------|---|
| Description | Recombinant Human Mucin-1 is produced by our Mammalian expression system and the target gene encodin... |
| Endotoxins | Less than 0.1 EU/ug |
| Form/Appearance | Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4 |
| Storage | Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months. |
| Tag | C-hFc |
| Note | For research use only |
| Uniprot ID | P15941 |
| Expiration Date | 6 months from date of receipt. |