

Product Datasheet

Human Lymphotoxin β R protein (orb737452)

Description	Recombinant Human Lymphotoxin beta Receptor is produced by our Mammalian expression system and the t...
Endotoxins	Less than 0.1 EU/ug
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-6His
Note	For research use only
Uniprot ID	P36941
Expiration Date	6 months from date of receipt.