

## Product Datasheet

### Mouse HVEM protein (orb737465)

<b>Description</b>	Recombinant Mouse Herpesvirus Entry Mediator is produced by our Mammalian expression system and the ...
<b>Endotoxins</b>	Less than 0.1 EU/ug
<b>Form/Appearance</b>	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
<b>Storage</b>	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
<b>Tag</b>	C-hFc
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">Q80WM9</a>
<b>Expiration Date</b>	6 months from date of receipt.