


Product Datasheet

UBE1 antibody (orb73754)

Description	UBE1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P22314
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:NT2D1 whole c...



Confocal immunofluorescence analysis of ...