



Product Datasheet

ERCC8(CSA) antibody (orb73763)

Description	ERCC8(CSA) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q13216
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:MOLT4 whole c...



Confocal immunofluorescence analysis of ...