



Product Datasheet

RFC4 antibody (orb73770)


Description	RFC4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P35249
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:H1299, B:HeLa



C:HepG2 whole cell lysate using RFC4 ant...



Immunohistochemical analysis of formalin...