



Product Datasheet

TOK alpha antibody (orb73799)


Description	TOK alpha antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9P287
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:A431, B:Hela ...



Merged with DNA probe



Confocal immunofluorescence analysis of ...