



Product Datasheet

Coronin 1B antibody (orb73807)

Description	Coronin 1B antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9BR76
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:293T whole ce...



Immunohistochemical analysis of formalin...