

Product Datasheet

Human Lipocalin-2/NGAL ELISA Kit (orb738505)

| | |
|----------------------------|---|
| Description | Human Lipocalin-2/NGAL ELISA Kit |
| Reactivity | Human |
| Range | 156 pg/ml - 10,000 pg/ml |
| Tested Applications | ELISA |
| Immunogen | Expression system for standard: NS0; Immunogen sequence: Q21-G198 |
| Storage | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| Note | For research use only |
| Sample Types | culture supernatants, serum, plasma (heparin), saliva and urine |
| Uniprot ID | P80188 |
| Sensitivity | 10 pg/ml |
| Expiration Date | Please enquire. |

