



Product Datasheet

CDC45L antibody (orb73867)

| | |
|----------------------------|---|
| Description | CDC45L antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | IF, WB |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| Note | For research use only |
| Clonality | Recombinant |
| Uniprot ID | O75419 |
| Expiration Date | 12 months from date of receipt. |



Western Blot analysis of A:H1299 whole c...



Confocal immunofluorescence analysis of ...