



Product Datasheet

Sec61 alpha antibody (orb73926)


Description	Sec61 alpha antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P61619
Expiration Date	12 months from date of receipt.



Confocal immunofluorescence analysis of ...



Western Blot analysis of A:K562 whole ce...



Immunohistochemical analysis of formalin...