



Product Datasheet

myotubularin related protein 4 antibody (orb74039)

| | |
|----------------------------|---|
| Description | myotubularin related protein 4 antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | IF, WB |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| Note | For research use only |
| Clonality | Recombinant |
| Uniprot ID | Q9NYA4 |
| Expiration Date | 12 months from date of receipt. |



Western Blot analysis of A:HeLa S3 whole...



Confocal immunofluorescence analysis of ...