



Product Datasheet

Cytokeratin 13 antibody (orb74040)


Description	Cytokeratin 13 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P13646
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:Hela whole ce...



Immunohistochemical analysis of formalin...



Confocal immunofluorescence analysis of ...