



Product Datasheet

Carbonic anhydrase 9 antibody (orb74096)

Description	Carbonic anhydrase 9 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q16790
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:MOLT4, B:Raji...



Confocal immunofluorescence analysis of ...