



Product Datasheet

HLA-DRB3 antibody (orb74135)


Description	HLA-DRB3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P79483
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:A549 whole ce...



Costained with Hoechst 33342



Confocal immunofluorescence analysis of ...