



Product Datasheet

20S Proteasome alpha6 antibody (orb74180)


Description	20S Proteasome alpha6 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P60900
Expiration Date	12 months from date of receipt.




Western Blot analysis of A:293T, B:A431,...



Immunohistochemical analysis of formalin...



Merged with DNA probe



Confocal immunofluorescence analysis of ...