



Product Datasheet

TFIID (SI) antibody (orb74193)

| | |
|------------------------|---|
| Description | TFIID (SI) antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| Note | For research use only |
| Clonality | Recombinant |
| Uniprot ID | P20226 |
| Expiration Date | 12 months from date of receipt. |



Western Blot analysis of A:293T whole ce...



Immunohistochemical analysis of formalin...