



Product Datasheet

succinate-CoA ligase, ADP-forming, beta subunit precursor antibody (orb74212)


Description	succinate-CoA ligase, ADP-forming, beta subunit precursor antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9P2R7
Expiration Date	12 months from date of receipt.




Western Blot analysis of A:HepG2, B:Molt...



Immunohistochemical analysis of formalin...



Costained with Hoechst 33342



Confocal immunofluorescence analysis of ...