

Product Datasheet

zona pellucida binding protein isoform 1 antibody (orb74233)

Description zona pellucida binding protein isoform 1 antibody

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications IHC, WB


Storage Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Note For research use only


Clonality Recombinant

Uniprot ID [Q9BS86](#)

Expiration Date 12 months from date of receipt.



Western Blot analysis of A:A431, B:HeLaS...



Immunohistochemical analysis of formalin...