



Product Datasheet

transmembrane emp24 protein transport domain containing 9 precursor antibody (orb74253)

Description	transmembrane emp24 protein transport domain containing 9 precursor antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9BVK6
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:H1299 whole c...



Immunohistochemical analysis of formalin...