



Product Datasheet

RNA binding motif protein 28 antibody (orb74258)

| | |
|----------------------------|---|
| Description | RNA binding motif protein 28 antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | IF, WB |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| Note | For research use only |
| Clonality | Recombinant |
| Uniprot ID | Q9NW13 |
| Expiration Date | 12 months from date of receipt. |



Western Blot analysis of A:293T whole ce...



Confocal immunofluorescence analysis of ...