



Product Datasheet

RAP2C, member of RAS oncogene family precursor antibody (orb74266)

Description	RAP2C, member of RAS oncogene family precursor antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9Y3L5
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:MOLT4 whole c...



Immunohistochemical analysis of formalin...