



Product Datasheet

C1orf165 antibody (orb74270)

Description	C1orf165 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q7L4P6
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:Molt-4 whole ...



Confocal immunofluorescence analysis of ...