



Product Datasheet

Yip1 domain family, member 4 antibody (orb74278)

| | |
|----------------------------|---|
| Description | Yip1 domain family, member 4 antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | IHC, WB |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| Note | For research use only |
| Clonality | Recombinant |
| Uniprot ID | Q9BSR8 |
| Expiration Date | 12 months from date of receipt. |



Western Blot analysis of A:A431 whole ce...



Immunohistochemical analysis of formalin...