



Product Datasheet

fibronectin type III and SPRY domain containing 2 antibody (orb74287)

| | |
|----------------------------|---|
| Description | fibronectin type III and SPRY domain containing 2 antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | IHC, WB |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| Note | For research use only |
| Clonality | Recombinant |
| Uniprot ID | A1L4K1 |
| Expiration Date | 12 months from date of receipt. |



Immunohistochemical analysis of formalin...



Western Blot analysis of A:H1299, B:Hela...