



Product Datasheet

LOC727787 antibody (orb74306)

Description	LOC727787 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P43630
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:A431, B:H1299...



Confocal immunofluorescence analysis of ...