



Product Datasheet

IRAK Antibody (orb74307)

Description	IRAK Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IP, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P51617
Expiration Date	12 months from date of receipt.



Western blot analysis of THP-1 (THP) and...



Immunohistochemical analysis of formalin...