


Product Datasheet

FLIP Antibody (orb74333)

Description	FLIP Antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	O15519
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH/3T3 whole c...



Immunohistochemical analysis of formalin...