



Product Datasheet

CX3CR1 Antibody (orb74362)

Description	CX3CR1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P49238
Expiration Date	12 months from date of receipt.



Western blot analysis of human spleen ly...



Immunohistochemical analysis of formalin...