



Product Datasheet

DcR3 Antibody (orb74377)

Description	DcR3 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	095407
Expiration Date	12 months from date of receipt.



Western blot analysis of human heart (H)...



Immunohistochemical analysis of formalin...