



## Product Datasheet

### IKK epsilon Antibody (orb74427)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | IKK epsilon Antibody  |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | ELISA, IHC, WB  |
| <b>Storage</b>             | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| <b>Note</b>                | For research use only   |
| <b>Clonality</b>           | Recombinant   |
| <b>Uniprot ID</b>          | <a href="#">Q14164</a>  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |



Western blot analysis of Jurkat whole ce...



Immunohistochemical analysis of formalin...