


Product Datasheet


IL-21 Antibody (orb74449)

Description	IL-21 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9HBE4
Expiration Date	12 months from date of receipt.



Immunocytochemical analysis of formalin-

...



Western blot analysis of HL-60 cell lysa...