



Product Datasheet

KappaB ras1 Antibody (orb74457)

Description	KappaB ras1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9NYS0
Expiration Date	12 months from date of receipt.



Western blot analysis of RAW264.7 cell l...



Immunohistochemical analysis of formalin...