



Product Datasheet

TLR2 Antibody (orb74498)

Description	TLR2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	O60603
Expiration Date	12 months from date of receipt.



Western blot analysis of A20 cell lysate...



Immunohistochemical analysis of formalin...