



Product Datasheet

TLR4 Antibody (orb74500)

Description	TLR4 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	000206
Expiration Date	12 months from date of receipt.



Western blot analysis of PC-3 cell lysat...



Immunohistochemical analysis of formalin...