



Product Datasheet

ACE2 Antibody (orb74520)

Description	ACE2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9BYF1
Expiration Date	12 months from date of receipt.



Western blot analysis of human kidney ly...



Immunohistochemical analysis of formalin...