

Product Datasheet

ATM (phospho-Ser1981) antibody (FITC) (orb7453)

Description	ATM (phospho-Ser1981) antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Tested Applications	FC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human ATM around the phosphorylation site of Ser1981
Target	phospho-ATM (Ser1981)
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q62388
Dilution Range	Flow-Cyt=1µg /test,IF=1:100-500
Expiration Date	12 months from date of receipt.