



Product Datasheet

TRIM5 alpha Antibody (orb74531)

Description	TRIM5 alpha Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q0PF16
Expiration Date	12 months from date of receipt.



Western blot analysis of human stomach (...)



Confocal immunofluorescence analysis of ...