



Product Datasheet

TLR11 Antibody (orb74540)

Description	TLR11 Antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q6R5P0
Expiration Date	12 months from date of receipt.



Western blot analysis of RAW264.7 cell l...



Immunohistochemical analysis of formalin...