



Product Datasheet

TAB1 Antibody (orb74568)

Description	TAB1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q15750
Expiration Date	12 months from date of receipt.



Western blot analysis of 3T3 cell lysate...



Immunohistochemical analysis of formalin...