

Product Datasheet

West Nile Virus Envelope Antibody (orb74589)

| | |
|----------------------------|---|
| Description | West Nile Virus Envelope Antibody |
| Species/Host | Rabbit |
| Reactivity | Virus |
| Conjugation | Unconjugated |
| Tested Applications | ELISA |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| Note | For research use only |
| Clonality | Recombinant |
| Uniprot ID | P06935 |
| Expiration Date | 12 months from date of receipt. |