

Product Datasheet

PS210 (orb746259)

Description	PS 210 is a potent, specific activator of PDK1 that target the PIF-binding pocket of PDK1, does not ...
Note	For research use only
Formula	$C_{19}H_{15}F_3O_5$
MW	380.319
CAS Number	[1221962-86-2]
SMILES	<chem>O=C(O)C(C(C1=CC=CC=C1)CC(C2=CC=C(C(F)(F)F)C=C2)=O)C(O)=O</chem>
Expiration Date	12 months from date of receipt.