

Product Datasheet

ML298 (orb746284)

Description	ML298 is a potent, selective phospholipase PLD2 inhibitor with IC50 of 355 nM, displays > 53-fold se...
Note	For research use only
Formula	$C_{22}H_{23}F_3N_4O_2$
MW	432.4388354
CAS Number	[1426916-02-0]
SMILES	<chem>FC1C=C(C=CC=1F)C(NCCN1CCC2(N(C3C=CC=C(F)C=3)CNC2=O)CC1)=O</chem>
Expiration Date	12 months from date of receipt.