

Product Datasheet

TTK21(CBP-p300 activator TTK21) (orb746291)

Description	TTK21 (CBP-p300 activator TTK21) is a small molecule activator of CBP/p300 histone acetyltransferase...
Note	For research use only
Formula	$C_{17}H_{15}ClF_3NO_2$
MW	357.7547143
CAS Number	[709676-56-2]
SMILES	<chem>C1C=CC(NC(=O)C2C=CC=CC=2OCCC)=CC=1C(F)(F)F</chem>
Expiration Date	12 months from date of receipt.